

New & Refurbished

# Veterinary Equipment

Supplies, Services & Support in Canada



*Medna*   
Medical Equipment Inc.

**Office BC: 604 283 4717**

**Office ON: 647 948 4757**

**Toll: 833 333 2581**

[sales@mednamedical.com](mailto:sales@mednamedical.com)  
[andrew@mednamedical.com](mailto:andrew@mednamedical.com)

# Top 4 Reasons Veterinarians Across Canada Choose Medna Medical Equipment.

1. The “Medna Customer Bill of Rights” ([See Page 32 Now for Full Details!](#))
2. Competitive Pricing, Expert Advice, Friendly Ongoing After Sale Support.
3. Medna is the Canadian Service Depot for 100% of the Products Represented  
(Avoid the Unnecessary, Frustrating, Service, Repair/Warranty Runaround).
4. No Cross-Border Nightmares... No Unexpected Customs Delays; No Costly Hidden Brokerage Fees  
([Ask us for a cost breakdown comparison direct from USA](#)).

Medna Represents The Top  
Veterinary Equipment  
Manufacturers & Suppliers  
Across Canada

---

## More Important Reasons Leading Veterinarians Across Canada Choose Medna

- Same Day Shipping even large capital items like X-Ray Dental Units, etc.,
- You start earning income from what you just paid for fast.
- Local Support - Offices in BC, Ontario Alberta and Ottawa.

---

## Did You Know Medna Specializes In Everything X-Ray?

- New Clinic Design, Specifications, Supply, Installation, Training, Support.
- HARP, RPS Surveys, Calibration, Repair.
- Digital Flat Panels, Service Repair Plus Maintenance Service Agreements.
- Canadian Compliant, Secure DICOM Cloud Backup System & Support.

For fast Local Support, Expert Equipment Advice, Local Service or Maintenance  
**Call Toll Free 833 333 2581 or Request Quote Online at [MENAMEDICAL.COM](http://MENAMEDICAL.COM)**

# MEDNA VETERINARY Contents

2	X-ray Systems
3	X-ray Solutions DR/CR
4,5	CT
6	C-Arm
7	X-ray Certificate
8	X-ray Accessories
9	Dental X-ray
10	Dental Software
11	Dental Care
12	Oxygen Concentrators
13,14	Lights
15	Monitors
16,17	Services
18	Ultrasound & ECG
19	Pumps
20	Autoclaves
21	Centrifuges
22	Warmers
23	Instruments
24	Supplies & Accessories
25	Consumables
26,27	Run Doors, Cages & Kennels
28,29,30,31	Tables: Exam & Surgery
32	Bill of Rights
33	Special Offers!

# Veterinary X-Ray System

Software  
Includes TTA  
VHS, TPLO  
Measurements

View your images anywhere, anytime,  
on any internet enabled device

- . Web PACS: plus 6 months unlimited cloud Storage and user access included.
- . Remotely View, store, send and receive your DICOM images 24/7/365.
- . Use a clinic iPad or tablet education & diagnosis review.

## 4-Way Float Top Table

- . Dimension:  
64 x 28"
- . Optional:  
SMART collimator with ATD (Automatic Thickness Detection)  
Easy clean-up surface with urine trap

Optional Instant  
Preview Capability



## Flat Panel Detector (5 Year Warranty)

- . 14x17 & 17x17 tethered or wireless
- . DIRECT Deposit (Csl AED) for low dose, high-quality diagnostic images, every time

## Easy Touch Intuitive Control Panel

- . Preset Factors: Or use the one-tap customization
- . Optional Image Preview: Stay table side and save time



# X-Ray Components

High - Voltage Cable



Collimator



X-Ray Tube



High Frequency  
Generator




\* These items can be sold separately or as a complete unit

# Digital X-Ray Solutions



## DR Flat Panel & CR Scanner

- . Superior image quality captured in 2-3 seconds
- . Viewing & Measurement Tools: Zoom, contrast, brightness, side-by-side, etc.
- . Measurements: TTA, VHS, TPLO
- . Full DICOM compliance with PACS software
- . Capable of transmitting images via e-mail and IP
- . Archival long-term storage and retrieval capabilities
- . On site & remote digital application training
- . Instant delivery
- . Exemplary technical support in Canada



## Upgrade to Digital on Any Existing X-Ray System Regardless of its Age In Only Half a Day!

### OPTIONS INCLUDE

- . TETHERD
- . WIPELESS
- . 10x12
- . 14x17
- . 17x17

### Low Dose Direct Deposit CsI, Second Generation AED with 5 Years Warranty



mini-PC



Equine Package



CR Scanner

- . Standalone, Industrial-grade computer
- . HIPPA - Compliant encryption
- . ARM - based security module
- . Medical grade work stations and monitors available

Go Digital For Only \$ 1.00  
No Payments up to 6 Months!

## Software for Image Acquisition, Processing, Diagnostics & PACS Storage

- . CaptureX imaging software with diagnostics.
- . OR technology for optimized workflow.
- . 20/20 Imaging Opal RAD
- . Nexus Varex - Nexus DR
- . DIRA Veterinary

# NewTom 5G XL

## ADVANCED 3D VETERINARY IMAGING

### Outstanding quality and innovation, exceptional features.

NewTom 5G XL, a CBCT device incorporating unique technology, guarantees advanced top-level diagnostics by maximising patient stability and providing a broad selection of FOVs, resulting in 3D and 2D images of the very highest quality.

CBCT technology combines low X-ray doses with high spatial resolution: the outstanding quality of 5G XL images allows veterinary specialists to identify pathologies and fractures that would otherwise be unidentifiable. Patient health - enhanced by the ECO Dose mode and exclusive SafeBeam™ technology - is a key design priority.


The native 21 x 19 cm FOV and numerous smaller high resolution FOVs constitute a toolbox that is perfect for veterinary diagnostics, which demand sharp and extremely detailed images.



Ask Us  
a bout our  
Refurbished  
Inventory



**SUPERIOR  
DIAGNOSTIC  
QUALITY**



**MOTOR-DRIVEN  
PATIENT TABLE**



**MINIMUM TIMES,  
MAXIMUM  
USER-FRIENDLINESS**



**VETERINARY  
SOFTWARE**


## SUPERIOR DIAGNOSTIC QUALITY

Superior standard  
3D examinations  
with a device  
designed for excellent  
performance.

Ultra-high definition volumetric images with native isotropic voxel resolution and minimal artifact presence. Thanks to CBCT technology, 5G XL offers fast scanning and extremely contained X-ray doses, increasing safety for both operator and patient and ensuring a more efficient workflow.

The outstanding versatility of 5G XL makes it the perfect device for otorhinolaryngological and muscular-skeletal diagnoses, with or without radiocontrast agents, for complete upper airway analysis and for the investigation of dental and maxillofacial pathologies.

### COMBINING TECHNOLOGY, PERFORMANCE AND SAFETY



5G XL offers outstanding results with low radiation doses for the patient, all thanks to undeniable constructive excellence:

The X-ray generator, the only one of its kind in the world, with rotating anode and a 0.3 mm focal spot, optimises performance by adapting emissions to the specific needs of the examination. It also provides high filtration levels to protect the animal from more harmful radiation types.

The large Flat Panel sensor with high signal-noise ratio improves final image quality, while collimation limits exposure to the specific regions of interest.


Innovative volumetric reconstruction algorithms allow the user to control the image chain and maximise diagnostic potential while minimising the presence of artifacts.

## MOTOR-DRIVEN PATIENT TABLE.

The NewTom 5G XL is characterised by a motor-driven, carbon fibre patient table. Controlled via a panel on the machine, it ensures comfortable, easy access to the gantry for operator and patient alike.



Motor-driven table control is even easier when performed via the PC as it allows the user to adjust patient positioning directly from the workstation. The open gantry makes access to the scan area easy and simplifies, along with the three reference lasers, the management of initial patient positioning (dorsal, ventral, rostral, cranio-caudal or caudo-cranial).



360° reconstruction




eXtra FOV vision



HiRes analysis

# Skanc C VET

## High Frequency Veterinary Mobile Surgical C-Arm



Ask Us

Mobile, High-Frequency Surgical C-Arm with 2.2kW Monoblock, 9" Tri-Mode Image Intensifier, kV Range of 40-110kV, 10.4" Touch-Screen Console with Image Adjustment Features, On-Board Storage For Up to 100,000 Images, 1k x 1k High Resolution Camera, Road Mapping & DSA (Standard), Laser Aimer, DICOM 3.0 with PACs integration, Dual 18" High Contrast Monitors and Sleek On Board Keyboard.

### Radiation:

- \*Safest C-arm in its class with near zero leakage radiation
- \*Lower dosage with excellent image quality

Generator:

- \*High-performance generator with up to 250kHz switching frequency

- \*Customized and in-house developed X-ray tube

### Modes:

- \*Pulse mode for a lower dosage
  - \*DSA & Roadmap comes as a standard feature
- Service support:

- \*Best-in-the-industry warranty terms
- \*Remote diagnostics capability

Low dosage patient positioning:

- \*Virtual collimation for ROI selection w/o X-ray
- \*Laser Aimer for initial patient positioning

Workstation capability:


- \*More than 100k images storage
- \*DICOM 3 compatible; PACS integration
- \*2-hour inbuilt power back-up
- \*Upgradable option for wifi capability for PACS
- \*Upgradable option for external monitor



# + Web PACS CaptureX Unlimited Viewers

## Reasons Why Veterinary Clinics Love Using The Medna Web PACS Solution:

- . Easy to Use
- . Automatic Backups!
- . Saves You & Your Staff Time and Money
- . Access your Studies in Seconds - From Anywhere!
- . Unlimited Viewing Stations
- . Works on Any Device (PC, Mac, iPhone, iPad, Android, Tablets etc...)



## X-Ray Certificate

### Mandatory Collage Requirement

- . Radiation Protection Survey (RPS) & HARP
- . Radiation Safety Certificate
- . Certificate of Compliance



**Go to Mednamedical.com**  
and book your appointment to  
receive 10% off

\* Please note: Other names are used for this test. Please contact us for details.

# Veterinary Portable & Equine X-Ray

## PXP20HF Plus & VET-20BT

- . 40-100 KV @ 1.0 kv increment
- . Lightweight
- . Small size
- . Freedom of use
- . Zero clutter
- . Point and shoot



Lightest, smartest and most durable portable X-Ray



# X-Ray Accessories

Leaded window frames



Lead Sheeting



Lead glasses



Lead Apron & Glove Rack



Lead Thyroid Collar



Lead Apron with Integrated Thyroid Collar



Veterinary Lead Apparel and Accessories Pack



Lead Mittens



Lead Gloves (One Pair)



# Dental X-Ray

Affordable & High-Quality DC Dental X-Ray Systems

In stock,  
next day  
shipping  
Call Today  
833 333 2581

## DCX5000W Dental X-Ray (Wall Mount)

Introducing a power pack and intelligent dental system

- . Low dosage
- . Sharp and crisp image quality X-Ray
- . Triggers from distance with no-shift image position
- . Enhancing doctor's experience and flexibility
- . Exclusive LCD display
- . Clean power withdrawal while on stand by
- . Energy conservation with shorter exposure times



Sensor

## DCX5000M Dental X-Ray (Mobile)

Portable version for safe, quick transport and convenient to use

Single unit with Multi features:

- . 0.4 - 0.5 mm focal point
- . Arm Reach: Short 60", Standard 72", Long 80"
- . Options: In-operatory or remote timer
- . Fully digital-programmable operation
- . Fast exposure system
- . Instant switching between small and large animals or film and digital sensor setting



Mobile Unit



Wall Mount


## Redefining Intraoral Sensor

- . Specifically designed for veterinary use
- . Fiber optic plate with over 2m cable length
- . Waterproof, can be immersed in disinfectant liquid
- . High resolution + High Sensitivity = Accurate Diagnosis
- . Powerful imaging software, 3sec image transfer
- . Enhanced image quality, fully DICOM, Email, Export

Buy a Wall Mount or Mobile Unit with the Digital Intraoral Sensor and

**SAVE \$500**

Why  
Choose  
DC?



DC units produce a more consistent x-ray output which allows for shorter exposure times and take quality images consistently, compared to AC machines.

# Dental Software



## captureX Dental Software

### Software Key Features

- . Developed in Canada
- . Cloud software
- . Intuitive and smart UI (User Interface)
- . Unlimited Storage



## Vault & View

- . Developed in Canada
- . 365/24/7 Cloud backup
- . Patient filters & Search
- . Wide variety of editing tools
- . Online access from anywhere



### Vault & View

To Register Visit:

[Mednamedical.com](http://Mednamedical.com)

# Oxygen Concentrators



## Nuvo Lite / 5

- Best for offices with only one anesthesia machine
- Years of trouble free operation
- Can be mounted on some gas machine carts
- Lockable control flow valve

<b>Oxygen Concentration</b>	90% (+6.5%/ -3%) at 5 LPM
<b>Liter Flow:</b>	0.125 to 5 LPM
<b>Sound Level:</b>	45 dBA typical
<b>Outlet Pressure:</b>	7.0 PSI
<b>Weight:</b>	32lbs/ 14.5kg
<b>Dimensions:</b>	14"Wx9" Dx23" H

**Available Accessories:**

- Mounting Bracket



## Nuvo 8/10

- Best for large practices
- Can accommodate up to 3 Anesthesia Machines
- Can be used with Oxygen therapy systems
- Oxygen supply for the whole clinic

<b>Oxygen Concentration</b>	90% (+5%/ -3%) at 8 LPM
<b>Liter Flow:</b>	2 to 8 LPM
<b>Sound Level:</b>	48 dBA typical
<b>Outlet Pressure:</b>	15 PSIG
<b>Weight:</b>	55.5lbs/ 25.2kg
<b>Dimensions:</b>	16"Wx16" Dx28" H

**Available Accessories:**

- Mounting Bracket

*Little  
Maintenance*

*Quiet  
&  
Reliable*

*Call  
Today!*

# Dental Care

## Scale-Aire

Always  
in  
Stock

Canadian  
Serviced

Support  
and  
Training

- High speed, multi-function veterinary dental unit
- All-in-one unit lets you easily switch from polishing to scaling to cutting.
- Surgical stainless steel, heat treated tips glide smoothly across tooth surfaces.
- Optional, fiber optics system with high-speed hand piece.
- Tips gently cool and lubricate teeth, and wash away dislodged debris
- Built-in 25K ultrasonic scaler maintains optimum frequency.



## Accessories & Maintenance Kits

Looking for a dental unit Maintenance Kit? Look no further we have maintenance kits for your Sonus V, Son-Mate II, Scale-Aire & Scale-Aire Mini's.

- Choose from: Basic - Deluxe packages.
- Spray Lubricants, One-Step conditioners, Drop Oil, Water Filtration Bottles, Cleaning Brushes, Prophy Angels & Bur Kits.
- Not seeing what you need? Contact us we will find it for you.



### Dental Scaler

- Water line with filter and quick disconnect.
- Perform curettage & deep curettage.
- Scale heavy calculus and remove stains.
- Transducer (stack) inserted in the hand piece.



### Scale -Aire Mini

- A powerful and potent tool against periodontal disease.
- Connects to a nitrogen tank easily.
- Water flows internally through the tip directly to the operative site, gently cooling and lubricating the tooth, washing away the dislodged debris.



### Son-Mate II

- From 200 RPM high torque micro-motor with forward/reverse switch.
- One foot control for both scaling/polishing modes.
- Includes a scaler hand piece with an inserted stack (transducer), assorted tips and a tip wrench.

# Veterinary Lighting

## AIM LED

- 45,000 lux (4181 fc) at 1 meter.
- Uses 55% less energy compared to the traditional light sources.
- 3 LED modules with an 20,000-hour life.
- Adjustable focus using center handles.



AIM - 100 and AIM - 50 Models Also available

### Philips Burton Guide to Lights for Veterinary Applications

Dept./Specialty	Good	Better	Best
Major Surgery	AIM-100	AIM-200 OR	APEX
Diag	Outpatient LED	CoolSpot II	AIM LED
Examination	Super Exam 50 LED	NOVA LED	Outpatient LED
Dental	Super Exam 50 LED	Gleamer / NOVA LED	CoolSpot II LED



LED Rated Life 50,000 hours

## Outpatient III LED

- 90,000 lux (7990fc) at 24" (61cm).
- The unique shape of the light head reduces shadows on the examination area
- 20,000 hours average LED life

#1 Dental & Treatment Choice



All Lights Available In Ceiling, Wall, Floor & Double Head



Dual Head



Ceiling Mount



Wall Mount



Mobile

### UV Light Magnifiers



- 2.25x total power magnification
- 8" focal length
- 6000 hour average bulb life
- 5- year limited warranty

### Nova Exam LED



- Low level: 25,000 lux (2323fc) at 18" (46cm).
- High level: 50,000 lux (4645fc) at 18" (46cm).
- Patented K-arm with a maximum arm reach of 46" (114.3cm) allows for quick and reliable positioning

### CoolSpot II



- Up to 60,000 lux (5574 fc) at 1 meter.
- Bi-level dimmer for variable intensity.
- Variable spot size of 2" - 10" (5cm-25cm).
- New, improved and easy to clean housing design with convenient Sterihandle for positioning

### Super Exam 50 LED



- 39,000 lux (3623fc) at 24" (61cm).
- LED light module with 40,000 hour life.
- Delivers value with longer life and lower energy consumption.
- Drift-free K arm with 42" (107cm) arm range ensures easy positioning.

## Medical illumination

### MI-300 MI-150

- MI 300 : Toggle between 30,000 and 15,000 LUX
- MI 150 : 15,000 LUX
- 3 Year Warranty
- MI 150, 300, 550, 750, 1000
- Ceiling, Wall OR Portable mounts
- 15,000 to 100,000 LUX models
- Outstanding quality
- Affordable price
- 50,000 Hours





# Monitors



## EPM 12V VET



## EPM 10V VET



### Reliable Quality

The ePM 12M Vet comes equipped with a high-quality shell materials, fan-free design and veterinary specific accessories. They are safe, durable, and built to last.

#### High-quality shell materials

- Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners.

#### Fan-free design

- Keeps fur out of equipment
- Reduces cross-contamination, protecting the health of patients and staff
- No internal dust, extending the life of the monitor

#### Veterinary specific accessories

- Low-trauma ECG clips maintain consistent contact
- Appropriately-sized NIBP cuffs provide accurate readings
- Durability-tested for long-lasting reliability

### Reliable Quality

The ePM 10 Vet comes equipped with a high-quality shell materials, fan-free design and veterinary specific accessories. They are safe, durable, and built to last.

#### High-quality shell materials

- Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners.

#### Fan-free design

- Keeps fur out of equipment
- Reduces cross-contamination, protecting the health of patients and staff
- No internal dust, extending the life of the monitor

#### Veterinary specific accessories

- Low-trauma ECG clips maintain consistent contact
- Appropriately-sized NIBP cuffs provide accurate readings
- Durability-tested for long-lasting reliability

# Accessories

## ECG CLIP



ECG: Non-traumatic crocodile clips

## NIBP CUFFS



NIBP: Variety of cuffs fits multiple patient sizes

*Medna*

Medical Equipment Inc.

## Services

For More Details & Request a Quote:  
Go to [MEDNAMEDICAL.COM](http://MEDNAMEDICAL.COM)



### Equipment Service

- Warranty (Medna is The Canadian Service Depot for 100% of products sold)
- Service Agreements
- Maintenance (with or without service agreement)
- Laser Equipment Service & Repair in Canada
- Loaner Equipment
- Exchange/Upgrade Programs



### X-Ray Service

- Radiation Safety Survey & Certificate RPS
- Calibration
- Repairs
- Digital X-Ray Support
- Cloud PACS Setup
- Lead Shielding & Power Requirements

**For Everything Veterinary In Canada!**

**Call 833 333 2581**

17

# MEDNA

is The Canadian Service  
Depot for 100% of products sold



## 7 Day-a-week Support in Canada

- Emergency Service 24/7  
(Next-day support guarantee)
- On-site & Remote Technical Support
- Repairs (Why replace it when you can fix it?)
- Retrofits/Upgrades  
(Take advantage of new technology)



## Financing

- Banking Relationships
- Leasing
- 3 or 6 months no payment plan

## Consulting

- New Clinic Design: Turn-key
- Room & Radiology Preparation
- Ergonomically friendly designs



## Training

- Staff Training
- Technique Optimization
- Digital X-Ray Training
- Dental X-Ray Training

## Installation

- New & Refurbished Equipment
- Retrofits & Upgrades

**For Everything Veterinary In Canada!**

**Call 833 333 2581**

# Ultrasound

## Vetus E7

Hand-Carried Veterinary Ultrasound System



Since its founding, Mindray Animal Medical has been dedicated to developing innovative and accessible ultrasound solutions for the animal care environment. It is this commitment to excellence that has brought you ZONE Sonography® Technology+ (ZST+), evolving industrial designs and comprehensive disinfection solutions. The Vetus E7 Ultrasound System has been thoughtfully designed to overcome the obstacles clinicians face in today's challenging healthcare environment. With a sealed user interface design and revolutionary, software-based beamformer technology, ZST+, the Vetus E7 System combines best-in-class image quality with an intuitive user experience to help ensure reliable and efficient diagnosis during the most challenging exams

## DP-50 Vet



With its streamlined features and compact design, the DP-50 Vet Ultrasound System is the ideal solution for a mobile veterinary environment. The System's powerful engine and high-speed processing provide outstanding image quality and a complete information management solution. For veterinarians seeking improved ergonomics and innovative ultrasound technologies, the DP-50 Vet is an excellent choice for providing high-quality care

## Z60 Vet



The Z60 Vet Ultrasound System is a comprehensive choice for veterinarians seeking innovative and high-performing ultrasound technology. A compact System, the Z60 Vet combines unsurpassed image quality with ease-of-use, and practical functionality, while meeting the financial goals of any practice.

## M6 Vet



### Intelligent Workflow

- iTouch™: Intelligent image optimization for B, color and PW mode
- Smart Tracking: Continuously track the color flow and optimize the best color box in real time scanning
- Smart Doppler: Automatic optimize the size and position of color box and doppler gate placement to ensure optimal signal color and spectral doppler sampling

### Comprehensive Functionalities

- iScape™: real-time panoramic imaging with velocity indicator, a very useful tool to evaluate the racing horse tendon
- Free Xros™ Imaging: Anatomical M mode, it can support EF measure offline and significantly shorten the exam time
- Advanced cardiac function assessment packages including Free Xros™ CM Imaging: Anatomical Curved M mode and

### Professional Design for Veterinary Use

- Veterinary information management tool for more than 6 types of animals
- Up to 22 veterinary body marks
- Abdomen, cardiac, OB, vascular, small parts and MSK exam modes and analysis packages for different animals
- Veterinary OB table
- User-defined formulas for more applications

## M9 Vet



### Dedicated Veterinary Software

- Data information
- Exam modes
- Bodymarks
- Measurements and reports

### Advanced Features

- High quality CW imaging
- TDI and quantitative analysis for confident cardiology analysis
- Tissue Tracking quantitative analysis for advanced myocardial study
- Auto measurement including Auto PW trace, Auto EF etc.

## DP-30 Vet



### Dedicated Veterinary Software

- Data information
- Exam modes
- Bodymarks
- Measurements and reports

### Advanced Features

- High quality CW imaging
- TDI and quantitative analysis for confident cardiology analysis
- Tissue Tracking quantitative analysis for advanced myocardial study
- Auto measurement including Auto PW trace, Auto EF etc.

# Pumps

## EN-V5 Infusion



- 4.3 inch color touch screen, parameter setting and editing directly
- Uses any standard IV sets and stores no less than 2000 drugs in library
- Multi-function operation with 5 infusion/inclusion modes
- Infusion rates range from 0.1 to 2000ml/h
- Last therapy for rapid infusion start
- 10 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP44 waterproof

## EN-V7P Infusion



- 4.3 inch color touch screen, parameter setting and editing directly
- Uses any standard IV sets and stores no less than 30 drugs in library
- Multi-function operation with 3 modes
- Electric door and anti-free flow clip patented design
- Double CPU for increased safety, Ultrasonic air-in-line detector
- Historic record more than 5000 logs
- 9 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP24 waterproof, and stackable design

## EN-S7P Syringe



- 4.3 inch color touch screen, parameter setting and editing directly
- Automatic recognition of syringe sizes (5, 10, 20, 30 and 50/60)
- Stores no less than 30 drugs in library
- Multi-function operation with 3 modes
- KVO range from 0.01 to 5ml/h, editable
- History record more than 5000 logs
- 12 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP24 waterproof, and stackable design

## VV-100 Vein Viewer

- Easy to operate
- Accurate imaging
- Three operating modes



**Basic Mode**  
Basic monochromatic mode,  
clear image



**Enhance Mode**  
Two-color mixing,  
enhanced appearance



**Inverse Mode**  
Inverse color image,  
promoting vein



## Tubing set

- Optimally deliver IV fluids, TPN solutions, blood and blood products with a +/- 3% accuracy.
- All control-V IV pumps are now calibrated to this tubing in channel 1.
- 103", 20 drops/ml, dual injection ports, rotating luer lock.
- 50 IV tubing sets/cases.

## IV Stands



# Autoclaves

## Clave23+ with Compressor Class S Autoclave, Non-Vacuum



Full automatic autoclave with LCD displaying, independent steam generator, Independent temperature and pressure monitoring, higher precision control System, automatic alarm for water level, door locking and error code Indicated. \*Includes Built in Printer and Datalogger  
# of Trays: 4 trays (standard equipped), max 5 trays

<b>Chamber Size</b>	F9.72"X17.7" (23L)
<b>Overall Size</b>	18.9"(W) x 18.3"(H) x 27.6"(D) (23L)
<b>Nominal power</b>	1600W
<b>Voltage</b>	110VAC - 60Hz

## Tuttnauer Manual Autoclaves



- Automatic shut off at the end of both sterilization and dry cycles.
- A long life electro polished Stainless Steel chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.

	10M	9M
<b>Chamber Size</b>	10"x18"	9"x18"
<b>Volume</b>	6gal/23L	5gal/19L
<b>Cassette Capacity</b>	3full & 3half	2full & 2half
<b>No. of Trays</b>	4	3



## Ultrasonic Cleaning Systems

- High-efficiency 60kHz cavitation for a thorough and quieter cleaning.
- Durable vinyl clad steel housing. Easy to clean and scratch resistant.
- The 1-gallon system (CSU1) holds a standard selection of instruments or 2 half sizes cassettes. The gallon holds 3 full size cassettes.

## Sterilization Ideals



Instrument Pouch Racks



Chamber Brite Autoclave Cleaner



Ultrasonic Cleaning Tablets



Water Distiller

# Centrifuges



## E8 Digital

The E8 Fixed is designed for blood separations at a fixed speed of 3500 rpm only (1534 g-force). The E8 Digital will separate blood at 3500 rpm and can be slowed to separate other fluids at lower g-force, such as urine specimens as outlined by C.L.I.A. regulations.



Best  
Seller

## E8 Combination

The E8 Combination Centrifuge is specifically designed to run 75mm microhematocrit tubes and 1ml-15ml micro tubes and test tubes. This centrifuge features a quiet operation for a more peaceful laboratory environment. ing speed function for precision results.



## E8 Portafuge

The portable E8 Portafuge is designed for spinning on the go. This model is loaded with features and improvements, such as higher speeds for faster and cleaner separations, quieter operation, lower chamber temperatures for cooler test-tube samples, and auto-calibrating function for precision speeds.

We carry a full line of Centrifuge, for more models please contact us!

# Microscopes



## Revelation III DIN

- LED illumination, for daylight color
- Low heat, low power consumption
- Super-long bulb life (30,000)



## i4 Infinity

- Exceptional optical quality
- Rotating head to adjust the height
- Cool, long-life LED illumination



## Innovation with BioVIEW

- Live presentations and client education
- Premium BioVID HD1080p up to 4K
- 1080p 11.6" HD Monitor or external 4K TV



MiniVID WiFi Microscope Camera



BioVIEW Camera and Monitor



BioVID 4K 8MP Ultra HD  
Microscope Camera

# Warmers



## Wild Cat Warmer

Compact unit delivers powerful, accurate warming while taking up minimal workspace in your facility. It has a large crisp display & is lightweight which is ideal for transport

- Less than 10 pounds
- Temperature sensing hose-end helps to deliver accurate air temperature every time.
- Four standard one-touch temperature selections: Ambient, 32°, 38°, 43°



## Bedjet Warmer

A set of 3 different reusable, machine washable warming blankets are included making it by far the best value on the market.

- Remote control
- Adjust from anywhere
- Weighs 13.8lbs
- Super Quiet Operation



## HTP-1500 Warmer

The HTP-1500 is the safest and most accurate localized heat therapy system available. Its easy-to-read display provides continuous real-time temperature readings. The digital controller uses proprietary software to manage your prescribed set-point temperature to within one-degree Celsius (1.0°C)!\* Three temperature safety limits guard against the possibility of overheating and the Hi Limit Switches can be tested without taking the unit apart! The HTP-1500 is also Mercury Free and uses Tap Water!

# Blankets & Pads



# Aspirators

## Gomco Optivac

### AC/DC Portable Suction Unit

- 25-550mm Hg vacuum range, coupled with greater than 30lpm flow.
- Rechargeable battery with 3 hour run time at maximum vacuum.



## Schuco - Vac

### Dependable Suction Unit

- Aspirator wire basket accommodates a wide variety of canisters
- Quick-connect tubing for all connections





# Instruments

## High-Quality, German Surgical Instruments

### Dental Instruments

- Surgical Blades
- Periotomes
- Elevators
- Root Tip Pick
- Extraction Forceps
- Needle Holders

### Surgical Packs

- Feline & Canine Spay packs
- Standard Veterinary Eye Pack
- Standard Veterinary Orthopedic Pack
- Small Animal Dental Extraction Kits
- Rabbit Dental Extraction Kits


 Medna  
Medical Equipment Inc.

## To Serve You Even Better... We now sell consumables

Because using high-quality, long-lasting and intelligently designed instruments make completing your procedures more efficient, safe and less stressful. We can provide your clinic or hospital with precision surgical instruments for every procedure at extremely competitive pricing.

Call with any important questions or to arrange a visit so you can see the instruments first hand. We distribute over 10,000 products consisting of, High Quality Surgical Implements & consumables.

#### Balfour Retractor


• 90mm • 180mm • 300mm



#### Baby Gelpi Retractor



#### Feline | Canine Spay Pack



# Supplies & Accessories

**Refurbished items available:**  
Please contact us  
for more details!

Stethoscopes



Wall-Mounted Diagnostic Sets



Welch Allyn Diagnostic Sets



Thermometer



Thermometer



Tono-Vet



Digital Touch screen Electro-Son



Apnea Alert



Frameless Stretcher



Cat Grabber



Stretcher



Positioners



Digital Platform Scale



Cat Scale



# Consumables

Sterile Gloves



Surgical Masks



Sundries



Loupes



Hydrogen Peroxide



CaviWipes / Solution



Scalpel Blades & Handles



Sutures



CarboLime



Syringes & Tips



Aftercare Cones



Suture Removal Kits



Laceration Kit



Puppy Pads



Sterilization Indicating Tape



Surgical Gowns



Surgical Drapes



# Run Doors & Cage Accessories

**Grill Run Door**



**Glass Run Door**



**Full Frame Run Doors with Grill**



**Full Frame Run Doors with Glass**



**Oxygen Therapy Door**



**Intensive Care Unit**



**Transport Cage with Therapy Door**



**Noise Reduction Strapping**



**Clipboards**



**Cage With Pole**



**Clipboards with Bin**



**Humidifier – Single Patient Use**



**Plastisol Removable Cage Floors**



# Cages & Kennels

All Styles of Regal Cages and Designed with these Quality Features:

**Doors** Stainless steel door spokes are deburred smooth to eliminate sharp ends. Stainless steel, electropolished removable ambidextrous door may be remounted with hinges on the right side. Door extends below cage floor to help keep paper liners in place. Door fits snug to the cage body to eliminate areas where an animal's paws could be pinched. 1.875-in. (.48cm) diameter stainless steel vertical spokes are spaced 1-in. apart and welded at every intersection to crossmembers and frame. 3.125-in. (.79cm) diameter stainless steel door frame and crossmembers are welded together for strength.

**Hinges** 3.125-in. (.79cm) diameter stainless steel hinge pins allow easy removal of door. Heavy 12-gauge stainless steel deformed barrel hinges operate smoothly for the life of the cage.

**Hardware** Stainless Steel installation hardware is the strongest available. Vertical molding strips and 14-gauge bolt-on corner plates mount easily with .250-20 stainless steel truss-head machine screws.



**Latch** Stainless Steel Latch easily disassembles for sanitizing and adjusting. .3125-in. (.79cm) diameter stainless steel latch keeper latches door securely. Heavy 12-gauge stainless steel latch brackets. Latch connector bar is notched for easy opening. Heavy 12-gauge stainless steel self-latching door latch connector is designed to accommodate a padlock for greater security.

**Body** 22-gauge Type 304 stainless steel body is heli-arc welded into one piece to create the strongest cage available. .750 in. (1.90cm) coved rear corners on top and bottom eliminate accumulation of dirt and bacteria. .750-in. (1.90cm) radius bends on all sides and back provides a smooth, easy-to-clean interior surface. Distinctive floor design with fluid-retaining front lip is standard on all Regal Cages. Bottom lip overlaps lower cage when used in a Regal Cage arrangement.

**Frame** 18-gauge stainless steel die-formed floor support brace keeps floor rigid for animal comfort. 18-gauge stainless steel die-formed rear support keeps cage level. Hidden 18-gauge stainless steel die-formed front frame makes the front of the cage extra strong.

## Sample Configurations



Single Door Stainless Steel Cages



Double Door Stainless Steel Cages



Aristo-Cat Kitty Suites, Cat Condo



Stationary



Heated



Mobile

# Exam Tables



**Classic Pedestal Exam Table with Tilt Top**



**Classic Pedestal Exam Table with Regal 300 Electronic Scale**



**Classic Pedestal Exam Table**

## Lift Table



*More models & configurations available. Contact us for details.*

## Premier Exam Table



## Classic Vet Mate Gurney




## Classic Tubular Base Exam Table



## Classic Wall-Mounted Exam Table



## Folding Wall-Mounted Exam Table



# Prep-Procedure, Wet & Dental Table

Regal Stainless Steel  
Imperial Laminated



Dental Model



Middle Space



Full Cabinet



Classic  
Model

Classic  
Model II



Classic  
Model I



V-Rack Top



Solid Top



Options  
Available

Multi Depth Sink  
Coated Removable  
Bathy Rod




Recessed Sink



Full Sink  
Stainless Steel Removable  
Bathing Rack



L-Arrangement,  
2station



# Surgery Operating Tables



Classic  
Pedestal



Elite  
Hydraulic



Regal  
Electric

## Our Selection Include:

Flat Top or V-Style Top Heated or  
Non-Heated 46" or 58" Wide Top

## 4 Easy Steps To Choosing Your Table:



V-Top or Flat-Top?



46" or 58"?



Heated or Non-Heated?



Hydraulic, Manual or Electric?



# Tables & Surgery Accessories

Over Instrument Table



Utility Tables



Over Operating Table



IV-Stand



Chrome Base Stool



Mayo Stands



5-Caster Instrument Stand



Solution Stand



Each Pull-out Tray is sold separately



Vet Stretcher



Sponge Receptacle



Kick Bucket



Bucket



Stainless Steel IV Poles



Stainless Steel Instrument Tables



Operating Table Casters



Head Extension for V-Top Table



Dual-Tilt Retrofit Kit for Classic & Elite Operating Tables



Thoracic Positioner



# Customer's Bill of Rights

**Medna Medical is proud abide by the Medna Customer's Bill of Rights for the Veterinarian Industry which Protects our customers from any negative experience they might otherwise have from other veterinarian supply companies:**

- 1 We will treat our customers the way we would like ourselves to be treated, like family!
- 2 We guarantee our customers will receive honest and ethical advice, service and pricing.
- 3 The customer's best interest always take precedence in any intransaction with Medna , and the customer is always right.
- 4 We will clearly explain the pros and cons of all options presented so the customer can make a wise, informed decision.
- 5 We will deliver on 100% of our promises. If we can not follow through on a commitment, we will not promise what we can not deliver!
- 6 The customer will always receive prompt communication from their representative.
- 7 Clear expectations will always be set and agreed upon by the customer and representative, to avoid misunderstandings and unnecessary problems.
- 8 All quotes will clearly state all prices and fees up front-protecting you from surprise costs like shipping, installation, duties, brokerage fees, and price changes.
- 9 Service vehicles will be fully stocked with the necessary diagnostic and repair equipment ensuring the quickest repair and shortest service times possible.
- 10 We will protect our customers with the best warranties in our industry.
- 11 We will invest in the most advanced equipment in our industry to ensure our customers get best and fastest service possible.
- 12 Our technicians will have no less than 5 years-experience before working in the field for our customers.
- 13 We will always remember to quote on all items required to have you purchase work as you expect it to do.

# Special Offers

Go Digital For  
**Only \$1.00** No Payments  
 for **6 Months!**

## Convert Your Existing X-Ray to Digital



**SPECIAL  
 DISCOUNT**

### Are You...

- Waiting more than 2-3 seconds for diagnostic quality images?
- Tolerating (or embarrassed by) poor image quality?
- Retaking X-Rays?
- Exposing yourself to more radiation than necessary...?

### Then NOW is the best time to upgrade and save thousands!

Choose From Three Special Offers Plus Discounted Lease Options:

- 60 Month term, 10% Buyout
- 66 Month term , \$1.00 Buyout
- Skip 3/6 months and

## START MAKING MONEY NOW!

Ask us about our Loyalty Program!

Get 1 Year Of Unlimited Web PACS and Cloud Storage, Unlimited Viewing Stations AND A 5-Year

**“Bumper to Bumper Warranty”**

# How Does Your Veterinary Supplier Rate?

For your convenience, we've provided this important checklist with all the questions you'll want to ask before buying the products you need.

## Company Status

- Has Offices, Inventory & Staff (Not 1-Man Show or Drop Shipper)
- Stock 80% of Items in Canada (Has Same/Next Day Shipping)
- Single Vendor Capable (Full Product Line & Support, including X-Ray)

## Canadian Vs. USA Supply & Service

- Local Service & Support Technicians (No "Technician Fly in Costs and Delays)
- Canadian Warranty (No Ugly Cross-Border Nightmares!)
- Canadian Service Depot (for All Products Sold)

## Shipping & Customs

- No Unknown Fees for Shipping or Customs Brokerage
- Pick any shipping company you want (Some Force Expensive ones!)

## Service & Support

- Next Day Support Guarantee
- Maintenance & Service Agreements on 100% of Products Sold
- Loaner Equipment Available (Avoid Costly Downtime)

	Medna vet	Big International	Online	US Dealer
Company Status	✓ ✓ ✓	✓ ✗ ✓	✗ ✗ ✗	✓ ✗ ?
Canadian Vs. USA Supply & Service	✓ ✓ ✓	? ? ✗	✗ ✗ ✗	✗ ✗ ✗
Shipping & Customs	✓ ✓	✓ ✓	✗ ✗	✗ ✓
Service & Support	✓ ✓ ✓	? ✓ ?	✗ ✗ ✗	✗ ✗ ✗



Office BC: 604 283 4717  
 Office ON: 647 948 4757  
 Toll: 833 333 2581  
 sales@mednamedical.com  
 Andrew@mednamedical.com  
 www.mednamedical.com  
 124 Garden Ave, North Vancouver, BC, V7P 3H2, Canada

