

T-CLASSIC 9

9-Liter Portable Autoclave

Fast. Compact. Portable. Simple. Cost Effective
making your process easy, quick and affordable.

The T-Classic 9 is a 9-liter lightweight, portable steam autoclave designed as a cost-effective sterilization solution. Fully automatic with an easy one touch operation, delivering short cycles, starting at **18 Minutes!***

Ideally designed to meet portability needs of small to medium sized practitioners or clinics.

The T-Classic 9 can also be used in combination with another autoclave or as a back up.

T-Classic 9 allows you to have your instruments sterilized and ready quickly!

*Cycle times may vary, depending on load and weight.



Tuttnauer performance with the simplicity and affordability your small practice needs



Fast

The T-Classic 9 offers competitively quick unwrapped cycle times, in as short as 18 minutes*. The fast point-of-use sterilization makes your turnaround a priority, increasing efficiency, while eliminating bottlenecks.

*Cycle times may vary, depending on load & weight.



Intuitive & Simple

The T-Classic 9 is a simple one step operated device. With one touch, the T-Classic 9 will begin an automatic cycle until completion. An electronic detector will turn off the power if water level is insufficient. The perfect autoclave for limited space.



Compact & Portable

The T-Classic 9 is the lightest tabletop autoclave available, weighing only 11.7 lbs. allowing for easy portability, while requiring minimal counter space.



Safe

A thermal fuse is in place to protect against overheating and electrical faults.



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel

Technical Specifications

T-Classic 9

Max Width / Height (overall) 13.4" / 13"

Max Net Weight 5.3 kg/11.7 lbs.

Chamber Diameter/Height 8.2"/9"

Chamber capacity 9 Liter / 2.4 Gallon

Max instrument length 9"

Max load weight 3kg/6.61 lbs.

Sterilizing temperature 250°F

Min sterilizing time 18 min.

Voltage 120V

Watts 1200W

Frequency 50-60Hz

Designed and manufactured to comply with the following standards and guidelines: ASME VIII Division 1, Medical Device Directive MDD 93/42/EEC, FDA 510(K) - Number K962903, CSA, TSSA Canadian Design Registration

More products:



Ultrasonic Cleaners



Water Distiller



T-EDGE Autoclaves

LTR330-0121/T-ClassicUSA/REV-1/10.22

Tuttnauer
Innovation · Legacy · Partnership

Medna
Medical Equipment Inc.

T: 833 333 2581 P: 604-283-4717 C: 647 948 4797 C: 778 881 4902

E: Info@mednamedical.com