

Valueklave™ 1730

# Valueklave™ 1730



## Speed, Quality & Value

- Fast cycle time, 11 minutes from a hot start (Standard Unwrapped Cycle Time)
- Efficient drying cycle for wrapped instruments
- Contemporary design with easy to read panel
- Economically priced a compact size

# Valueklave™ 1730

The Valueklave™ 1730 will satisfy all of your basic sterilization needs at an economical price. High speed and reliability are built into every Valueklave™ 1730.



| Valueklave™ 1730 Specification      |   |
|-------------------------------------|---|
| Chamber Dimensions<br>0 x L in (mm) | 7"x13"<br>(170 X 340)                   |
| Volume                              | 2 gal/7.5 L                             |
| Overall Dimensions<br>DxWxH in (mm) | 17.9"x 17.4 "x 12"<br>(455 x 440 x 305) |
| Standard Cassettes                  | 2 half                                  |
| Tray Dimensions<br>DxWxH in (mm)    | 11.6"x 4.7 "x 0.8"<br>(295 x 120 x 20)  |
| No. of Trays                        | 3                                       |
| Standard Unwrapped<br>Cycle Time    | Cold: 16 min<br>Hot: 11 min             |
| Voltage (V)<br>Freq. (Hz)           | 120<br>50/60                            |
| Power (W)<br>Curret (A)             | 1350<br>11                              |
| Shipping Weight<br>lbs (kg)         | 55 (25)                                 |

## Features:

- Fast Cycle Time. 11 minutes from a hot start\* (Standard Unwrapped Cycle time).
- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Overheat cutoff prevents burnout of heating elements in case of insufficient water in the chamber.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.



\*Come up time, sterilization exposure and exhaust



Clean&Simple™ Ultrasonic Cleaners and Tablets



Steam Distiller



Chamber Brite Autoclave Cleaner

LT18330-0005 /VALF/USA/V-3 /10.16



T: 833 333 2581 P: 604-283-4717 C: 647 948 4797 C: 778 881 4902

E: [Info@mednamedical.com](mailto:Info@mednamedical.com)